

# Application Note:

Scale-up of *Staphylococcus aureus* bacteriophage production in the CellMaker Regular single-use bioreactor system.

QUALITY | RELIABILITY | SERVICE EXCELLENCE | PRODUCTIVITY

## Introduction


There is much commercial interest in exploiting the use of bacteriophages to control bacterial infections (phage therapy), with clear potential for control of wound infection and also general sterilization applications.

Novolytics, a company whose aim is to develop phage therapies to treat a number of problematic bacterial infections, is initially targeting Meticillin-Resistant *Staphylococcus aureus* (MRSA). Currently shaker flask cultures are used to progress these research programs, but clearly as progress is made this shaker flask scales will become limiting. Novolytics therefore evaluated the Cellexus CellMaker Regular (8L) to assess simple scale-up of *Staphylococcus aureus* bacteriophage production.

The protocol developed for this application dictates the need to seed the CellexusBag with *S. aureus*, then allow the bacteria to progress to early exponential growth phase before inoculating with bacteriophage. The phage-enriched media is then later harvested following successful phage infection and cell lysis.

## Experimental

The CellexusBag was filled with 8L LB media (+ 2mls SE-15 antifoam) and pre-warmed in the CellMaker enclosure before inoculation with 0.5% (40mls) *S. aureus* (overnight culture). Upon attaining an OD600 of 0.3, bacteriophage was introduced into the CellexusBag at a multiplicity of infection (MOI) of 0.25. The phage-enriched media was later harvested following phage infection and cell lysis (associated with a drop in OD600). These experiments are run at relatively low ODs and therefore the oxygen demand of the culture is not high. Typical run parameters are shown in the table (below).



TIME (mins)	OD600	Set Temp	Process Temp	Gas Mix (Air/O <sub>2</sub> /Headspace) (all L/min)	Pressure (mBar)	Comments
30	0.097	37	37	5/0.5/0.4	40	
70	0.162	37	37	5/0.5/0.4	40	
90	0.269	37	37	5/0.5/0.4	40	
95	-	37	37	5/0.5/0.4	40	1% phage (80ml) added
100	0.350	37	37	5/0.5/0.4	40	OD taken straight after phage addition
185	0.358	37	37	5/0.5/0.4	40	Sign of cell lysis
255	0.046	37	37	5/0.5/0.4	40	Culture cleared; harvest



Increasing protein expression yields of a 44kDa His-tagged protein in *Pichia pastoris* in a single-use bioreactor (CellMaker Regular) rapid, proof-of-concept experiment.

QUALITY | RELIABILITY | SERVICE EXCELLENCE | PRODUCTIVITY

## Introduction

*Pichia pastoris* is a strain of yeast that is very often preferred to *Saccharomyces cerevisiae*, other yeast strains, and *E. coli* because it generates more human-like glycosylation patterns/structures within recombinant proteins. As a result it has become a major method of recombinant human protein expression. Expression is driven by the tightly regulated AOX1 promoter, which is induced by methanol. Secretion of recombinant protein (directed by a signal sequence) into the media simplifies the downstream protein purification process.

The CellMaker Regular is a single-use bioreactor offering considerable ease-of-use and is designed to give high protein yield in oxygen-hungry bacterial and yeast expression systems. The system is available in both 8 litre and 50 litre formats. The following results were generated using the smaller system.

## Experimental details

The target protein (44kDa) contains a C-terminal His-tag. It is secreted into the supernatant using alpha-factor signal sequence. The expression levels are very low (<50 ug/L) when produced in shaker flasks. A modified BMMY media was used and supplemented by the addition of 1% MeOH and 0.1% glucose every 24 hours.

The table (below) outlines the settings and readings from the Cellexus Cellmaker during this experiment:


TIME (hours)	OD600 [Shaker]	Set Temp	Process Temp	Gas Mix (Air/O <sub>2</sub> /Headspace)	Pressure (mBar)
0	0.94 [0.94]	29	29	4/4/1	40
24	40 [33]	29	29	2/4/1	50
48	49 [33.8]	29	29	2/4/1	50
72	50 [32.5]	29	29	2/4/1	50
96	54 [33]	29	29	2/4/1	50

## Results

The graph below shows the increased OD and relates to the yield of biomass compared to shaker flasks prior to harvesting from a 6 litre fermentation.

In this sample enhanced aeration was achieved using an oxygen enriched gas supply (from the built-in oxygen generator) at 4l/min oxygen and 2l/min of air. Higher ODs resulted in 1.5-fold increase in biomass and 2-fold increase in protein yield (as determined by western blot analysis). Protein was found to be glycosylated correctly and was secreted appropriately into the medium prior to purification.

*Pichia pastoris* growth at 29 °C



## Conclusions

The data generated here shows that the 2-fold increase in expression levels of protein can be quickly achieved by switching from shaker flasks to a single-use bioreactor in a proof-of-concept of this type. With additional optimisation using this bioreactor it is likely that the yield of protein will be further enhanced. Scaling from the 8 litre unit to the identically designed 50 litre unit may offer a further 6-fold increase in yield.

*Data kindly provided by Dr. Maxim Rossmann, Laboratory of Molecular Biology, Cambridge, UK.*

For further details, or to request a quotation, contact us now.

## Contact details

Headquarters - United Kingdom  
 Cellexus International Limited  
 6 Riverside Court, Mayo Avenue  
 Dundee, DD2 1XD  
 United Kingdom

Email: [sales@cellexus.com](mailto:sales@cellexus.com)  
 Office: +44 (0)1382 666357

For full details of our distributors please see our website.

[www.cellexus.com](http://www.cellexus.com)

